





Collection Kit		Sample Type	LAB9913 (CT/NG)	LAB9912 (CT)	LAB9928 (NG)	LAB9964 (TV)	LAB9957 (BV)	LAB9958 (CV/TV)	LAB9960 (MG)	LAB9961 (BV/CV/TV/MG)	LAB9934 (HSV 1&2)	Aptima HPV	Aptima HPV Genotype
ThinPrep Pap Vial		Cervicovaginal	X	X	X	X						X	X
		Anal/rectal											X
Aptima Multitest Swab		Vaginal (Clinician Collect)	X	X	X	X	X	X	X	X	X		
		Vaginal (Patient Collect)	X	X	X		X	X	X	X			
		Throat	X	X	X								
		Anal/rectal	X	X	X								
		Anogenital Lesion										X	
Aptima Unisex Swab		Endocervical	X	X	X	X			X				
		Male Urethral	X	X	X	X			X			X	
		Anogenital Lesion										X	
Aptima Urine Specimen		Female Urine	X	X	X	X			X		X		
		Male Urine	X	X	X	X			X		X		

Test Order Code	Directory Link	Test Name	Acceptable Collection				Sample Type
			ThinPrep	Aptima MultiTest	Aptima Unisex	Aptima Urine	
LAB10332* or Pap Order	https://apps.aurora.org/dos/public/lab/lab10332	ThinPrep Image-Guided Liquid Based Pap Test (GYN)	X				Cervicovaginal
		ThinPrep Image-Guided Liquid Based Pap Test with reflex to Aptima HPV					
		ThinPrep Image-Guided Liquid Based Pap Test and Aptima HPV					
		ThinPrep Image-Guided Liquid Based Pap Test and Aptima HPV with reflex to HPV Genotype 16, 18/45					
LAB9939 or HPVHR	https://apps.aurora.org/dos/public/lab/lab9939	Aptima HPV	X				Cervicovaginal, anal/rectal, vulvar
LAB9938 or HPVGT	https://apps.aurora.org/dos/public/lab/lab9938	Aptima HPV Genotype 16, 18/45	X				Cervicovaginal, anal/rectal, vulvar
LAB9053 or CYTOLOGY	https://apps.aurora.org/dos/public/lab/cytology	Cytology Non-Gynecological	X				Anal/rectal, vulvar

*Note Chlamydia, Gonorrhea and Trichomonas can be added

Test Order Code	Directory Link	Test Name	Acceptable Collection				Sample Type
			ThinPrep	Aptima MultiTest	Aptima Unisex	Aptima Urine	
LAB9913 or CGRNA	https://apps.aurora.org/dos/public/lab/lab9913	Chlamydia trachomatis / Neisseria gonorrhoeae by NAA	X	X	X	X	Vaginal (clinician or patient collect), Endocervical, Male Urethral, Throat, Rectal, Urine
LAB9912 or CTRNA	https://apps.aurora.org/dos/public/lab/lab9912	Chlamydia trachomatis by NAA	X	X	X	X	Vaginal (clinician or patient collect), Endocervical, Male Urethral, Throat, Rectal, Urine
LAB9928 or NGRNA	https://apps.aurora.org/dos/public/lab/lab9928	Neisseria gonorrhoeae by NAA	X	X	X	X	Vaginal (clinician or patient collect), Endocervical, Male Urethral, Throat, Rectal, Urine
LAB9964 or TVRNA	https://apps.aurora.org/dos/public/lab/lab9964	Trichomonas Vaginalis by NAA	X	X	X	X	Vaginal (clinician collect only), Endocervical, Male Urethral, Urine
LAB9960 or SWOMG	https://apps.aurora.org/dos/public/lab/swomg	SwabOne Mycoplasma genitalium by NAA		X	X	X	Vaginal (clinician or patient collect), Endocervical, Male Urethral, Urine
LAB9957 or SWOBV	https://apps.aurora.org/dos/public/lab/swobv	SwabOne Vaginosis by NAA (Lactobacillus species, Gardnerella vaginalis and Atopobium vaginae - 1 test result)		X			Vaginal (clinician or patient collect)
LAB9958 or SWOCN	https://apps.aurora.org/dos/public/lab/swocn	SwabOne Candida/Trichomonas Panel by NAA (Candida species, C. glabrata, Trichomoniasis - 3 test results)		X			Vaginal (clinician or patient collect)
LAB9961 or SWOPNL	https://apps.aurora.org/dos/public/lab/swopnl	SwabOne Vaginitis Panel by NAA (Bacterial vaginosis, Candida vaginitis, Trichomoniasis, Mycoplasma genitalium - 5 test results)		X			Vaginal (clinician or patient collect)
LAB9934 or HSVPCR	https://apps.aurora.org/dos/public/lab/lab9934	Herpes Simplex Virus 1 & 2 by NAA		X	X	X	Anogenital lesion, Vaginal (clinician collect only), Male Urethral, Urine